




Tracking Molecular Conformations of Stearic Acid in Surface-Adsorbed Monolayers


Elizabeth L. Stewart, Xixuan Guo, Flaminia Marrucci
Nathan Price, J. Clayton Baum, Joel A. Olson

Florida Institute of Technology, Department of Chemistry, 150 West University Blvd., Melbourne, FL 32901





INTRODUCTION

Stearic acid was observed to adsorb onto highly-oriented pyrolytic graphite (HOPG) in 2 different conformations.





Horizontal
Dark
Features

Vertical
Dark
Features



Syn
Conformer

Anti
Conformer




More negative electrostatic charge is present near the carbonyl, thus it appears darker in the image.

- Gas-phase DFT shows the syn conformer is favored by 13.5 kJ/mol.
- Would result in a ~99:1 ratio at thermal equilibrium @298 K.
- STM image above shows ~2:1 ratio.
- Suggests a large barrier to rotation.

The purpose of this investigation is to track conformations over time to see if interconversion can be observed.

EXPERIMENTAL



A drop of solution of stearic acid in 1-phenyloctane was placed onto the HOPG surface.

The STM tip is immersed directly in the stearic acid solution to image the monolayer at the solution-HOPG interface.

- STM was used to observe stearic acid adsorbed onto the HOPG surface.
- Images were collected under low drift ambient conditions.
- Images were collected in constant current mode with a sample bias of +0.80 V and a setpoint of 100 pA.
- Surface was imaged, and features (conformers) were tracked for ~2 hours.

RESULTS

Conformers were identified by the shape of the carboxyl ring in the STM image.

syn conformer
anti conformer
syn conformer





Image size: 7 nm by 7 nm
 $V_{\text{bias}} = +0.80 \text{ V}$; $I_{\text{tun}} = 100 \text{ pA}$

Individual conformers were then tracked over time to observe whether or not they interconverted from anti to syn or vice-versa.

The following images show examples of molecules tracked during time-lapse STM imaging.



Of eleven conformer groups tracked in the experiment, there were no clear changes of conformation.

CONCLUSIONS

- Molecular conformations do not frequently interconvert.
- This suggests a high energy barrier to rotation.
- Adsorption mechanism must allow for either conformer to be present upon formation of the monolayer.

ACKNOWLEDGMENTS



The authors gratefully acknowledge the National Science Foundation (grant #1058427) for their generous support.

This research was reported in:
X. Guo, F. Marriucci, N. Price, E.L. Stewart, J.C. Baum, J.A. Olson,
J. Phys Chem. C, **2015**, 119, 24804–24811.